


Twin Lock-BCP Postero Medial Proximal Tibial Plate, 3.5mm



Indications

Proximal Medial Tibia Fixation

Screw Size

Twin Lock-BCP Screw \varnothing 3.5mm
 Twin Lock-BCP Cancellous Screw \varnothing 3.5mm
 Cortical Bone Screw \varnothing 3.5mm
 Twin Lock-BCP Screw \varnothing 3.7mm Cannulated

Matching Instrument

Code 6297.00 Instrument set for 3.5mm Locking system.

*manufactured on demand with MOQ.

Size (Holes)	Product Reference Number							Length in mm
	Titanium (G5)	Titanium (G5) Rev.1	Titanium (G23)	Stainless Steel (SS316L)	Stainless Steel (SS316L) Rev.1	Stainless Steel (SS316LVM)	Stainless Steel (SS316LVM) Rev.1	
02 holes	*2286.02	*C0240.2002	C0240.4002	1286.02	C0240.1002	*3286.02	*C0240.3002	
04 holes	*2286.04	*C0240.2004	C0240.4004	1286.04	C0240.1004	*3286.04	*C0240.3004	
06 holes	*2286.06	*C0240.2006	C0240.4006	1286.06	C0240.1006	*3286.06	*C0240.3006	
08 holes	*2286.08	*C0240.2008	C0240.4008	1286.08	C0240.1008	*3286.08	*C0240.3008	
10 holes	*2286.10	*C0240.2010	C0240.4010	1286.10	C0240.1010	*3286.10	*C0240.3010	
12 holes	*2286.12	*C0240.2012	C0240.4012	1286.12	C0240.1012	*3286.12	*C0240.3012	